



DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



MARKED

F
U
N
C
T
I
O
N
S

- **No edge frame**, in order to enable trunk rotation and side-arm motion. It's adjustable in height by knobs, without using wrenches, "tools free".
- **Birch wood table**. It is removable and adjustable in height, tilt and depth.
- **Trunk Support** adjustable in its circumference thanks to Velcro straps, in height, back and forward.
- **Pelvic Support** adjustable in its circumference thanks to Velcro straps, in height, back and forward.
- **Numbered bars** to find the right setting, determined by the therapist.
- **Knee-pads** adjustable in height, width, depth and in a rotatory movement.
- **Heelrests** adjustable in depth and width.
- **Birch wood footboard**. It gives a safe guide to the child.
- **Twinned solid Wheels**, 75 mm. diam., rear wheels with brakes. Thanks to the wheels the aid can be moved from one room to another having the child on it.
- **Electronic Kit**: patented system fixed under the footboard which can transform an "audio" electronic signal into vibration.

Please check off the interested products:



size 1


size 2

size 3



Attention!
Read carefully the user and maintenance manual

M
E
A
S
U
R
E
S



User dimensions

USER HEIGHT (M)
 size 1: min cm 75 - max cm 100
 size 2: min cm 90 - max cm 125
 size 3: min cm 120 - max cm 150

IMPORTANT:
make sure that measures F and G indicated in the figure do not exceed the suggested range indicated for the supports.

size	Max. Load Cap.	Weight	A	B	C	D*	E	F	G	L
1	Kg. 35	Kg. 24,6	66 cm	76-88 cm	49-77 cm	64 cm	14-24 cm	41-55 cm	53-72 cm	60 cm
2	Kg. 45	Kg. 26	70 cm	78-90 cm	59-87 cm	65 cm	14-36 cm	52-68 cm	67-94 cm	65 cm
3	Kg. 55	Kg. 28,5	73 cm	80-92 cm	69-100 cm	70 cm	14-42 cm	62-86 cm	82-110 cm	70 cm

C
O
M
P
O
N
E
N
T
S

944 Circumference reduction parts Thickness: 3 cm. each side

view from above



865 Multiadjustable headrest



859 Tablet holder



Notes:-----

